

### ARTICLE OPEN



# Direct identification of *ALK* and *ROS1* fusions in non-small cell lung cancer from hematoxylin and eosin-stained slides using deep learning algorithms

Chen Mayer 1,4 Efrat Ofek<sup>1,4</sup>, Danielle Even Fridrich<sup>1</sup>, Yossef Molchanov<sup>1</sup>, Rinat Yacobi<sup>1</sup>, Inbal Gazy<sup>2</sup>, Ido Hayun<sup>2</sup>, Jonathan Zalach<sup>2</sup>, Nurit Paz-Yaacov<sup>2</sup> and Iris Barshack<sup>1,3</sup>

© The Author(s) 2022

Anaplastic lymphoma kinase (ALK) and ROS oncogene 1 (ROS1) gene fusions are well-established key players in non-small cell lung cancer (NSCLC). Although their frequency is relatively low, their detection is important for patient care and guides therapeutic decisions. The accepted methods used for their detection are immunohistochemistry (IHC) and fluorescence in situ hybridization (FISH) assay, as well as DNA and RNA-based sequencing methodologies. These assays are expensive, time-consuming, and require technical expertise and specialized equipment as well as biological specimens that are not always available. Here we present an alternative detection method using a computer vision deep learning approach. An advanced convolutional neural network (CNN) was used to generate classifier models to detect ALK and ROS1-fusions directly from scanned hematoxylin and eosin (H&E) whole slide images prepared from NSCLC tumors of patients. A two-step training approach was applied, with an initial unsupervised training step performed on a pan-cancer sample cohort followed by a semi-supervised fine-tuning step, which supported the development of a classifier with performances equal to those accepted for diagnostic tests. Validation of the ALK/ROS1 classifier on a cohort of 72 lung cancer cases who underwent ALK and ROS1-fusion testing at the pathology department at Sheba Medical Center displayed sensitivities of 100% for both genes (six ALK-positive and two ROS1-positive cases) and specificities of 100% and 98.6% respectively for ALK and ROS1, with only one false-positive result for ROS1-alteration. These results demonstrate the potential advantages that machine learning solutions may have in the molecular pathology domain, by allowing fast, standardized, accurate, and robust biomarker detection overcoming many limitations encountered when using current techniques. The integration of such novel solutions into the routine pathology workflow can support and improve the current clinical pipeline.

Modern Pathology (2022) 35:1882-1887; https://doi.org/10.1038/s41379-022-01141-4

#### INTRODUCTION

Lung cancer is the leading cause of cancer-related mortality and accounts for ~1.76 million deaths per year worldwide. Non-small cell lung cancer (NSCLC) comprises ~85% of all lung cancers and typically presents at advanced stages. The introduction of genotype-directed targeted therapies for NSCLC has transformed the care of patients and dramatically improved survival<sup>1</sup>.

Several tyrosine kinase inhibitor (TKI) targeted therapies are currently used for the treatment of lung cancer. The current College of American Pathologist (CAP)/International Association for the Study of Lung Cancer (IASLC)/Association for Molecular Pathology (AMP) guidelines recommend screening advanced-stage lung cancer patients for targetable alterations including testing for *EGFR*, *ALK*, and *ROS1*<sup>2</sup>. The National Cancer Center Network (NCCN) guidelines go even further and recommend broad molecular profiling using multiplex mutation screening assays or next-generation sequencing (NGS)<sup>3</sup>.

Despite rapid progress in the field, many practical issues often prevent or delay the initiation of targeted therapies. These obstacles include insufficient tissue for testing, poor quality of DNA or RNA, sequencing failure  $(3-10\%)^{4,5}$ , high cost, and high turnaround times<sup>6</sup>. Effective approaches to reducing these barriers are vitally needed to ensure patients' access to available treatments.

Anaplastic lymphoma kinase (ALK) belongs to the insulin-like tyrosine kinase (TK) receptor superfamily. Its rearrangements occur in ~5% of NSCLC cases and are associated with younger age, no or light smoking history, and characteristic histology such as signetring cells and cribriform pattern<sup>7</sup>. ROS proto-oncogene 1 (ROS1) gene belongs to the subfamily of TK insulin receptor genes. Its rearrangements occur in ~1-2% of patients with NSCLC and are also commonly found in younger patients, no or light smokers, and are associated with histological features of adenocarcinoma and similar characteristics to ALK-fusions such as signet-ring and cribriform pattern<sup>8,9</sup>. While there are features associated with these alterations, none are in high enough association nor specific enough to be used as markers for fusions. Thus, all NSCLC patients are currently tested using conventional laboratory tests such as IHC, FISH and NGS as discussed in the manuscript. Both ALK and ROS1 rearrangements result in fusion proteins with constitutive TK activity. A variety of genes have been identified as their fusion

<sup>1</sup>Institute of Pathology, Sheba Medical Center, Ramat Gan, Israel. <sup>2</sup>Imagene AI, Tel Aviv, Israel. <sup>3</sup>Sackler faculty of Medicine, Tel Aviv University, Tel Aviv, Israel. <sup>4</sup>These authors contributed equally: Chen Mayer, Efrat Ofek. <sup>⊠</sup>email: chen.mayer@sheba.health.gov.il

Received: 21 March 2022 Revised: 20 July 2022 Accepted: 20 July 2022

Published online: 3 September 2022

partners, most frequently *EML4* for *ALK* and *CD74* for *ROS1*<sup>10</sup>. *ALK* fusions can effectively be blocked by TKIs, including alectinib, lorlatinib, brigatinib, crizotinib, and ceritinib<sup>3</sup>. Since *ALK* and *ROS1* are closely related kinases, many of these TKIs are also used for the treatment of *ROS1*-positive NSCLC patients<sup>11</sup>.

The initial gold standard method for detecting ALK and ROS1 rearrangements was fluorescence in situ hybridization (FISH). Immunohistochemistry (IHC) was later added as an accepted alternative to FISH for ALK and as a screening method for ROS1, requiring subsequent confirmation of positive cases by another molecular or cytogenetic method<sup>12,13</sup>. Additional methods for identification of ALK and ROS1-rearrangements include reverse transcription PCR (RT-PCR) and NGS<sup>14</sup>. Each of these methods has its advantages and disadvantages, and the concordance between them varies between 70.7-98% depending on the gene and the assays used<sup>2,15</sup>. FISH, while still the gold standard, is expensive, timeconsuming, and requires technical expertise and specialized equipment. Moreover, both false-negative and -positive cases of ALK and ROS1- rearrangements have been documented 16,17. The use of IHC to detect ALK rearrangements has reported sensitivities and specificities ranging from 95% to 100%<sup>2</sup>. However, several studies have described discrepancies between FISH and IHC methods, with <90% concordance for ALK and only 80% concordance reported for ROS1 testing 15,18. The challenges facing both NGS and RT-PCR methodologies include a high technical failure rate (as much as 10–30% for NGS)<sup>19</sup>, partly due to the poor DNA/RNA quality in some formalin-fixed paraffin-embedded (FFPE) samples, as well as other limitations such as tissue scarcity and high cost, which hinder the widespread clinical implementation of these methods.

Deep learning (DL) is a subset of Artificial intelligence (Al) that uses machine learning (ML) algorithms in artificial neural networks to detect complicated patterns in large datasets. Computer vision is a field of AI that trains computers to interpret and understand the visual world using digital images. DL has revolutionized computer vision and has become a standard technique for image classification and is increasingly used to improve clinical practice across the medical imaging domain in a wide range of applications such as CT, MRI and mammography<sup>20</sup>. The advances in digitalization of tissue slides, alongside the evolving utilization of DL in the biomedical field, together with the increasing demand for molecular profiling in cancer, have led to the recognition of the potential for DL use in cancer classification and prediction of biomarker status. Oncogenic driver mutations are most probably reflected in an altered morphological pattern of the cancer cells and/or tissue. While not visible to the human eye, these changes can be detected by visionbased algorithms and are distinctive to each specific alteration. Indeed, there is an increasing number of publications describing the application of DL in the detection of molecular changes in histological slide images. Examples include, among others, prediction of PD-L1 (n = 82 test cases)<sup>21</sup>, EGFR, KRAS, and TP53 (n = 59 test cases)<sup>22</sup> alterations in NSCLC, MSI in colorectal cancer (n = 479 test cases)<sup>23</sup>, and ER/PR/HER2 status in breast cancer (n = 2611-2714test cases)<sup>24</sup>. While these studies demonstrate the potential of Al in the molecular diagnosis of cancer, the results are not yet robust enough for the technique to replace the currently used genetic and molecular tests. To prove the robustness of the method, further studies should be performed on larger and more diverse clinical cohorts in order to demonstrate comparable performance to the standard of care.

Here we present the use of an Al-based solution to detect genomic rearrangements in *ALK* and *ROS1* from automatically scanned hematoxylin and eosin (H&E) stained lung tumor slides. Our work highlights the advantages of this Al solution, which provides a rapid, cost-effective, and accurate testing alternative. Integration of such solutions within the clinical practice can support pathologist workflow and assist oncologists to tailor the best treatment for optimal management of patient care, bridging the current gap towards true implementation of precision medicine.

# MATERIALS AND METHODS Sample collection

Digitized images of H&E slides from advanced-stage NSCLC patients were collected from the Pathology Department of Sheba Medical Center (according to IRB 7451-20-SMC). All images were scanned at 40× magnification using the Philips IntelliSite Ultra-Fast scanner (Philips Digital Pathology Solutions, Best, Netherlands) and converted into TIFF format.

#### Data

A total of 234 NSCLC cases tested for *ALK* and *ROS1* rearrangements at Sheba Medical Center between 2012–2021 were selected for the study. These were randomly divided into training (N = 162) and validation (N = 72) sets. *ALK and ROS1* status was determined based on IHC, FISH, or NGS. Approximately 80% of patients in the validation set were tested by at least two testing methods (Supplementary Table S1). The training set included H&E scanned images as well as *ALK* and *ROS1* status annotations, while the validation set included images only.

#### Fixation and H&E staining

Tissues were processed according to standard pathological procedures. Tissues were fixed in buffered formalin, embedded in paraffin, and sectioned using a microtome. Sections were then placed on histologic slides and stained using H&E.

#### ALK (D5F3) and ROS1 immunostaining

For the immunohistochemical studies, 4µm wide sections were prepared from FFPE blocks and positive control was added at the edge of the slides. IHC staining was performed using *ALK* (clone D5F3, V790-4794, Ventana Medical Systems Inc., Oro Valley, AZ, Oro Valley, AZ, USA) and *ROS1* (1:50, 3287S, Cell Signaling, Danvers, MA, USA) antibodies with the OptiView Amplification Kit (Ventana Medical Systems Inc.) in conjunction with the OptiView detection kit (Ventana Medical Systems Inc.) on a Benchmark Ultra staining module (Ventana Medical Systems Inc.). All cases were tested by a single and dedicated thoracic pathologist. As indicated by the International Association for the Study of Lung Cancer (IASLC) intensity scoring was as follows: strong staining (3+) is clearly visible using 2× or 4× objective lens, moderate staining (2+) requires a 10× or 20× objective lens and weak (1+) cannot be seen until a 40× objective is used. Cases with a diffuse staining and (3+) score were interpreted as positive.

#### ALK and ROS1 molecular analysis

For *ALK* and *ROS1* fusion analysis, DNA and RNA were extracted from paraffinembedded sections using the KingFisher kit according to the manufacturer's instructions. The concentration of extracted DNA and RNA was measured using NanoDrop 2000c Spectrophotometer (Thermo Fisher Scientific, Waltham, MA, USA) following DNA measurement using a Qubit Fluorometer (Thermo Fisher Scientific) and brought to a concentration of 25 ng/µl.

cDNA was synthesized from RNA using the lon AmpliSeq<sup>TM</sup> HiFi Mix (Thermo Fisher Scientific). Sequencing was performed on the lon Torrent System (Thermo Fisher Scientific).

#### **FISH analysis**

*ALK* and *ROS1* rearrangement detection was performed using Vysis LSI *ALK* (2p23) Dual Color Break Apart Rearrangement FISH Probe (Abbott Molecular, Abbott Park, IL, USA) and Zytolight SPEC *ROS1* Dual color Breakpart Probe (ZytoVision GmbH, Bremerhaven, Germany). Assays were performed as described previously<sup>25</sup>.

#### Algorithm development and validation

To achieve high model accuracies, the development of the DL model was performed in two steps: (1) unsupervised learning training to produce initial weights for the neural network architecture; and (2) algorithm fine-tuning using semi-supervised learning. Following model generation, a retrospective blinded validation was performed (Fig. 1).

Unsupervised learning training. 21,299 H&E-stained, FFPE pan-cancer, untagged whole slide images (WSI) were used from Imagene-Al LTD.'s internal database. All slides were de-identified and scanned at 40× magnification (0.23–0.27 µm per pixel resolution). No additional data was linked to the images. Data augmentation was applied to images using image pre-processing algorithms and generative adversarial networks




Fig. 1 Al algorithm development and validation. A A two-step approach was used to generate an ALK/ROS1 classifier; first, initial weights were generated from WSIs from Imagene-Al's internal database, using an unsupervised learning algorithm. This was followed by a fine-tuning step performed using a semi-supervised algorithm, specifically for the tested biomarkers. B A retrospective set of de-identified scanned H&E images of 72 NSCLC samples were collected and analyzed by the ALK/ROS1 classifier. Results of the Al solution and the routine practice (IHC, NGS, and/or FISH) were compared for performance assessment. FFPE Formalin-Fixed Paraffin-Embedded, WSI whole image slide.

(GANs). For this step, convolutional neural network (CNN) training was performed with an average of 87,000 image tiles ( $1536 \times 1536$  pixel) per WSI extracted using an on-the-fly image tile selection method (implemented to achieve random sampling of regions without the need for a priori creation of patches) fed through the network.

Algorithm fine-tuning. 162 H&E WSI cases of NSCLC biopsies were retrieved from Sheba Medical Center's archives and scanned at 40× magnification (0.25 µm per pixel resolution) using a Philips Ultra-Fast Scanner (Philips Digital Pathology Solutions). The training set included 15 ALK- and 7 ROS1-positive cases identified by gold standard testing (IHC, FISH and/or NGS). Data augmentation was applied on images using image pre-processing algorithms and GANs. To create the ALK/ROS1 classifier algorithm, fine-tuning of the initial neural network algorithm developed in the previous step was conducted using semi-supervised learning and a five-fold cross-validation scheme. First, DL algorithms were trained to detect regions of interest (ROIs) in the WSIs (Supplementary Fig. S1). Next, model training was conducted on image tiles (384 × 384 pixels) extracted from the ROI. For image tile extraction an on-the-fly image tile sampling method was used, with an average of 1 M tiles per WSI generated and used for the training. In a few instances, when a couple of slides from the same cases were available, the first serial slide was used. Retrospectively, we tested all the consecutive slides, with 100% concordance.

Algorithm validation. For the evaluation of the ALK/ROS1 classifier algorithm accuracy, NSCLC biopsies processed at Sheba Medical Center were used. A cohort of 72 de-identified WSI retrieved from Sheba Medical Center's archives was used as a double-blinded set. The algorithm status predictions were compared to the ground truth diagnosis, and specificity and sensitivity were calculated to assess the Al-based model's performance.

#### Inference of molecular status

Images were uploaded into the Al solution powered by Imagene-Al LTD. in a local station as de-identified H&E scanned images, without additional clinical information.

A categorical prediction was made by the ALK/ROS1 classifier model.

## RESULTS

#### ALK and ROS1 detection by conventional methods

In this study, we compared the performance of *ALK* and *ROS1* conventional testing methods to that of the Al solution. IHC, FISH and NGS were used as gold standard methods for the analysis of these genes' rearrangements. In total, 53 of the 72 cases in our cohort were subjected to NGS testing, 14 underwent FISH testing for *ALK*, and 7 for *ROS1* (Supplementary Table S1). *ALK* IHC was performed on 70 cases and *ROS1* IHC on 65 cases. Six out of 72 cases were detected as *ALK*-positive and two out of 68 cases as *ROS1*-positive. Five of the six *ALK*-positive cases were detected as positive by FISH, with one confirmed by NGS. The sixth case was identified by a combination of IHC and NGS. The two *ROS1*-positive cases were detected by FISH analysis.

#### Validation of the ALK/ROS1 classifier

The ALK/ROS1 classifier algorithm was developed as described in the methods. Briefly, models were developed and fine-tuned on a training set which included H&E scanned images together with ALK and ROS1 status annotations only. To validate the model, a retrospective cohort of 72 samples was tested with the AI solution using de-identified images data only (a representative heatmap of a patch for each alteration is presented in Supplementary Fig. S2). We then compared the results of the algorithm predictions with the results of the conventional methods described above. The model correctly inferred the status of all cases except for a single ROS1 false positive. This sample was analyzed by all traditional methods (IHC, FISH, and NGS), which did not identify any ROS1 alteration. Altogether, the model achieved a sensitivity of 100% for both genes and specificity of 100% and 98.57% for ALK and ROS1 respectively (Table 1). The negative predictive value (NPV) for ALK and ROS1 was 1, and the positive predictive value (PPV) was 1 for ALK-fusion and 0.505 for ROS1-fusion.

Table 1. Summary of ALK/ROS1 classifier results.

		Conventional methods		Al-based model						
	N	# Positive	# Negative	TP	TN	FP	FN	Sensitivity	Specificity	Concordance
ALK	72	6	66	6	66	0	0	100%	100%	100%
ROS1	68	2	66	2	65	1	0	100%	98.48%	98.53%

TP true positive, TN true negative, FP false positive, FN false negative.

Traditionally, AI model performance is measured using AUC (Area Under the ROC Curve) values. However, here we chose to evaluate our model's performance using the same characteristics as those employed to evaluate the assays used in clinical practice. The clinical performance of such assays (e.g., sequencing, IHC and FISH) is described by their sensitivity, specificity, and concordance with the current gold-standard tests. We thus reported the results accordingly; however, in accordance with the high performances, the AUC values were high as well, 1 for *ALK*-fusion and 0.986 for *ROS1*-fusion.

Of note, approximately 7% of samples analyzed by IHC (N=9) were reported as "equivocal" and four samples tested by NGS failed (7.3%) due to technical issues such as poor DNA/RNA quality and insufficient diagnostic material. For diagnosis, these samples required additional testing with an alternative method. In contrast, all the tested samples in the cohort (including the aforementioned samples) were successfully analyzed by the AI pipeline.

#### Method comparison

To examine the performance of the traditional methods and the Al solution, we compared the outcomes of the various methods. When comparing the conventional methods, two samples had conflicting results, negative by IHC and positive by FISH analysis. These cases were given an official result of *ALK*-positive according to the FISH analysis, which was in concordance with the *ALK/ROS1* classifier calling.

Comparison of the AI solution with the routine methods identified only a single case with conflicting results. In this case, IHC, FISH, and NGS analysis did not detect any rearrangement, while the *ALK/ROS1* classifier reported this case as *ROS1*-positive.

#### **DISCUSSION**

In this study, we used an Al-based solution to detect ALK and ROS1 fusions in NSCLC cancer patients. To the best of our knowledge, this is the first study that has successfully applied an Al-solution on clinical samples to identify ALK and ROS1 rearrangements directly from H&E pathological WSI of lung cancer patients. Despite the low frequency of alterations in these oncogenes (5 and 1-2% of NSCLC cases for ALK and ROS1, respectively), their detection is of high importance, as it determines clinical treatment<sup>11</sup>. Both IHC and FISH were initially used to detect these rearrangements. However, advances in technologies in parallel with the continued increase in the number of actionable biomarkers, have led to the incorporation of NGS panels in clinical diagnostics. Unfortunately, all these assays have limitations, and testing rates are variable in the clinical practice, with cost, quality standards, and access being substantial barriers<sup>6,19</sup>. Furthermore, the current ASCO/CAP/IASLC/ AMP guidelines' recommendations refer to a turnaround time of ten working days between sample receipt and reporting of the molecular test<sup>26</sup>, yet in practice most institutes have a much longer turnaround time of several weeks from the pathology review to the final molecular report<sup>27</sup>

One of the major challenges in lung cancer screening is the limited amount of patient specimens. In 10–20% of cases tissue is insufficient for comprehensive sequencing or analysis by alternative methods<sup>28</sup>. This will most probably be less of a challenge to Al classifiers, as reflected in this study, where, unlike conventional methods, all the


samples in the cohort received a biomarker status calling. The median number of patches per slide that went through the Al analysis was 266 (ranged: 4–7055) per case. The ability to infer biomarker status from as little as four patches per slide is of major importance in this context and implies that Al-based solutions may be less sensitive to limiting amounts of biopsied tissues.

Routine molecular testing in lung cancer presents several challenges, some of which we faced in our cohort. The first, is the commonly encountered technical failure of NGS reactions. In our validation set, 5.5% of samples that underwent sequencing failed. Secondly, IHC screening can sometimes have ambiguous results. Approximately 7% of our validation set was "equivocal" in the initial IHC screening, requiring additional tests. Thirdly, the use of multiple tests (FISH, NGS, and IHC) may sometimes result in conflicting interpretations. Indeed, relatively high levels of discordance between tests have been reported<sup>2,15,18</sup>. In our validation set, six of the eight positive cases were tested by at least two methods. Of these, two cases had conflicting results (negative in IHC and positive in the FISH assay). Two additional cases were equivocal in IHC and positive according to the FISH analysis (one for ALK and one for ROS1). The remaining two cases showed concordance between methods, representing only a third of the positive cases. This issue is of major importance, especially in positive cases, where identification of ALK or ROS1 fusions dictates targeted treatment regimens.

Of the 140 tests performed using the Al solution (72 for ALK and 68 for ROS1), only a single case showed discordance with the final pathology report, demonstrating high accuracy of the Al solution. Further testing will be required to ascertain the robustness of this assay, and to see whether it will maintain its accuracy in real-world scenarios with larger datasets and when implemented in additional clinical settings and medical centers. Nonetheless, the diagnostic accuracy of the Al solution presented here, which was at least comparable if not even higher than that of the gold-standard tests, suggests that Al-based molecular diagnostic solutions may be appropriate as clinical diagnostic tools in the near future.

The technology used here examined the existence of *ALK* or *ROS1* fusions in NSCLC, irrespective of their partnering genes. This is in contrast to NGS for example, which can capture the identity of the fusion partner as well. Current targeted treatments for *ALK* and *ROS1* rearrangement disregard the partnering genes' identity and there is limited information regarding any difference in response between them<sup>7</sup>. Thus, we generated a model that detects *ALK* and *ROS1* fusions, regardless of their precise partner. However, each partner can result in a fusion product with slightly different expression and/or activity. Such differences might be identifiable by Al models if deemed to be required, but this will necessitate a significantly larger dataset.

Integration of an AI solution within routine clinical pathological workflow holds great advantages. Figure 2 illustrates the integration of the AI solution directly after the scanning of the H&E slide enabling the biomarker status calling within minutes. The suggested workflow can expedite the time that elapses between biopsy and commencing targeted therapy by assisting pathologists in flagging samples that potentially carry alterations in real-time to prioritize their examination. Moreover, after



## **Integration Into Clinical Practice**

**Fig. 2 Integration of Al-based solutions in the clinical pipeline.** Schematic overview of the integration of the Al-solution into the routine pathological biomarker diagnostic practice. After the biopsy is cut and stained with H&E, it is digitally scanned and transferred for pathological examination. Samples are further processed by the traditional route including IHC screening as well as FISH and NGS, to assess *ALK* and *ROS1* status. A final biomarker report is produced after an average of four weeks. The Al software (Imagene-Al. LTD) automatically analyses the WSI directly after the scanning of the H&E slide, taking an average of 2 min for *ALK/ROS1* status inference. Flagging of samples carrying *ALK* or *ROS1* alterations can be viewed immediately after this process by the pathologist and assist in their further analysis prioritization. After the results of the routine test panel are ready, they can be compared to the Al classifier inference. Discrepant cases can then be subjected to additional testing and clarification. Al artificial intelligence, WSI whole image slide, IHC Immunohistochemistry, FISH Fluorescence in situ hybridization, NGS Next-generation sequencing.

completion of the conventional testing, samples that have discordance in their biomarker status between the routine tests and the AI solution can undergo additional confirmational tests before reporting back to the patient. Hence, the AI solution can improve biomarker detection accuracy and assist or even change patients' testing prioritization.

Applying AI in clinical diagnosis can offer standardization and reproducibility, compared to current assays that suffer from subjectivity and wide variation between different labs and kits. Integration of AI medical software can also expand the availability of biomarker analysis in oncology as they do not require the expertise and specialized equipment that traditional assays demand, while preserving clinical standards of test accuracy. However, as opposed to molecular and biochemical assays that measure specific and known parameters, the outcome of the AI solution inference is based upon factors that are complex, obscure and undefinable. As technology progresses and the use in clinical practice is extended, new features will need to be developed to elucidate the morphological patterns used by the algorithm for analysis.

Advances in oncology and the gradual adoption of digitization of WSIs have paved the way for the incorporation of AI tools in digital pathology. This approach can improve our ability to tackle complex processes involved in cancer development. The concept of inferring molecular alterations directly from the morphometric phenotypes of tissue specimens presents a significant advantage compared to other analytical tools. This highlights a major advantage of AI-based methodology, such as the convolutional neural network (CNN) used here, which learns to identify patterns associated with the inferred alteration without the need for any prior annotations (such as histopathological features). By identifying different alterations that affect the cells and/or tissue morphology, DL algorithms can predict which tumors carry clinically significant changes without requiring a priori knowledge

regarding the nature of the alteration (*i.e.*, fusion partner, breakpoints, etc.). This can be applied to a wide range of biomarkers associated with different conditions such as mutations, protein expression, signatures, additional structural variants and more.

The molecular and IHC methods described in this study to diagnose ALK and ROS1 alterations are used in most countries for advanced cases only. Among other reasons, cost and the overload of molecular exams since their introduction to oncology contribute to the limited number of patient samples ultimately sent for testing. Those eligible for molecular testing often wait several, excruciating weeks for results, time that can be crucial in advanced cancer stages. An Al-based solution such as that validated in this study, is a fast application, requiring no additional tissue, and can be widely implemented potentially for all newly diagnosed NSCLC patients regardless of stage. Furthermore, automatization of such analysis can reduce potential human and technical errors. Moreover, utilizing a single histological slide image for assessment of multiple biomarkers can align with the increasing decline in bio-specimen availability for molecular testing. Broad implementation can be relatively easily applied and while this may pose new clinical dilemmas, it should ultimately result in benefit for the patients.

#### **DATA AVAILABILITY**

The datasets generated during the current study are not publicly available and are available from the corresponding author, CM, upon reasonable request.

#### REFERENCES

1. Blumenthal, GM, Karuri, SW, Zhang, H, Zhang, L, Khozin, S, Kazandjian, D, et al. Overall response rate, progression-free survival, and overall survival with targeted

- and standard therapies in advanced non-small-cell lung cancer: US Food and Drug Administration trial-level and patient-level analyses. *J Clin Oncol* **33**, 1008-1014 (2015).
- Lindeman, NI, Cagle, PT, Aisner, DL, Arcila, ME, Beasley, MB, Bernicker, EH, et al. Updated molecular testing guideline for the selection of lung cancer patients for treatment with targeted tyrosine kinase inhibitors: guideline from the College of American Pathologists, the International Association for the Study of Lung Cancer, and the Association for Molecular Pathology. Arch Pathol Lab Med 142, 321-346 (2018).
- Ettinger, DS, Wood, DE, Aggarwal, C, Aisner, DL, Akerley, W, Bauman, JR, et al. NCCN guidelines insights: non-small cell lung cancer, version 1.2020. J Natl Compr Canc Netw 17, 1464–1472 (2019).
- Steeghs, EMP, Groen, HJM, Schuuring, E, Aarts, MJ, Damhuis, RAM, Voorham, QJM, et al. Mutation-tailored treatment selection in non-small cell lung cancer patients in daily clinical practice. *Lung Cancer* 167, 87-97 (2022).
- Kunimasa, K, Matsumoto, S, Nishino, K, Nakamura, H, Kuhara, H, Tamiya, M, et al. Improvement strategies for successful next-generation sequencing analysis of lung cancer. Future Oncol 16, 1597-1606 (2020).
- Ahn, MJ. Molecular testing in lung cancer: still big gap in implementation for realworld use. J Thorac Oncol 15, 1399-1400 (2020).
- 7. Kim, H & Chung, J-H. Overview of clinicopathologic features of ALK-rearranged lung adenocarcinoma and current diagnostic testing for ALK rearrangement. Transl Lung Cancer Res 4, 149-155 (2015).
- 8. Bergethon, K, Shaw, AT, Ou, S-HI, Katayama, R, Lovly, CM, McDonald, NT, et al. ROS1 rearrangements define a unique molecular class of lung cancers. *J Clin Oncol* **30**, 863-870 (2012)
- 9. Pan, Y, Zhang, Y, Li, Y, Hu, H, Wang, L, Li, H, et al. ALK, ROS1 and RET fusions in 1139 lung adenocarcinomas: a comprehensive study of common and fusion pattern-specific clinicopathologic, histologic and cytologic features. *Lung Cancer* 84, 121-126 (2014).
- Remon, J, Pignataro, D, Novello, S & Passiglia, F. Current treatment and future challenges in ROS1- and ALK-rearranged advanced non-small cell lung cancer. Cancer Treat Rev 95, 102178 (2021).
- Dong, J, Li, B, Lin, D, Zhou, Q & Huang, D. Advances in targeted therapy and immunotherapy for non-small cell lung cancer based on accurate molecular typing. Front Pharm 10, 230 (2019).
- Planchard, D, Popat, S, Kerr, K, Novello, S, Smit, EF, Faivre-Finn, C, et al. Metastatic non-small cell lung cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up. *Ann Oncol* 29, iv192-iv237 (2018).
- 13. Kalemkerian, GP, Narula, N, Kennedy, EB, Biermann, WA, Donington, J, Leighl, NB, et al. Molecular testing guideline for the selection of patients with lung cancer for treatment with targeted tyrosine kinase inhibitors: American Society of Clinical Oncology Endorsement of the College of American Pathologists/International Association for the Study of Lung Cancer/Association for Molecular Pathology Clinical Practice Guideline Update. J Clin Oncol 36, 911-919 (2018).
- Conde, E, Rojo, F, Gómez, J, Enguita, AB, Abdulkader, I, González, A, et al. Molecular diagnosis in non-small-cell lung cancer: expert opinion on ALK and ROS1 testing. J Clin Pathol 75, 145-153 (2022).
- Reguart, N, Teixidó, C, Giménez-Capitán, A, Paré, L, Galván, P, Viteri, S, et al. Identification of ALK, ROS1, and RET fusions by a multiplexed mRNA-based assay in formalin-fixed, paraffin-embedded samples from advanced non-small-cell lung cancer patients. Clin Chem 63, 751-760 (2017).
- Jang, JS, Wang, X, Vedell, PT, Wen, J, Zhang, J, Ellison, DW, et al. Custom gene capture and next-generation sequencing to resolve discordant ALK status by FISH and IHC in lung adenocarcinoma. J Thorac Oncol 11, 1891-1900 (2016).
- Davies, KD, Le, AT, Sheren, J, Nijmeh, H, Gowan, K, Jones, KL, et al. Comparison of molecular testing modalities for detection of ROS1 rearrangements in a cohort of positive patient samples. *J Thorac Oncol* 13, 1474-1482 (2018).
- Mok, T, Peters, S, Camidge, DR, Noe, J, Gadgeel, S, Ou, SI, et al. Outcomes according to ALK status determined by central immunohistochemistry or fluorescence in situ hybridization in patients with ALK-positive NSCLC enrolled in the phase 3 ALEX study. J Thorac Oncol 16, 259-268 (2021).
- Smeltzer, MP, Wynes, MW, Lantuejoul, S, Soo, R, Ramalingam, SS, Varella-Garcia, M, et al. The International Association for the study of lung cancer global survey on molecular testing in lung cancer. J Thorac Oncol 15, 1434-1448 (2020).
- Hosny, A, Parmar, C, Quackenbush, J, Schwartz, LH & Aerts, H. Artificial intelligence in radiology. Nat Rev Cancer 18, 500-510 (2018).
- Sha, L, Osinski, BL, Ho, IY, Tan, TL, Willis, C, Weiss, H, et al. Multi-field-of-view deep learning model predicts nonsmall cell lung cancer programmed death-ligand 1 status from whole-slide hematoxylin and eosin images. J Pathol Inf 10, 24 (2019).
- Coudray, N, Ocampo, PS, Sakellaropoulos, T, Narula, N, Snuderl, M, Fenyö, D, et al. Classification and mutation prediction from non-small cell lung cancer histopathology images using deep learning. *Nat Med* 24, 1559-1567 (2018).

- Yamashita, R, Long, J, Longacre, T, Peng, L, Berry, G, Martin, B, et al. Deep learning model for the prediction of microsatellite instability in colorectal cancer: a diagnostic study. *Lancet Oncol* 22, 132-141 (2021).
- Rawat, RR, Ortega, I, Roy, P, Sha, F, Shibata, D, Ruderman, D, et al. Deep learned tissue "fingerprints" classify breast cancers by ER/PR/Her2 status from H&E images. Sci Rep. 10, 7275 (2020).
- McCoach, CE, Le, AT, Gowan, K, Jones, K, Schubert, L, Doak, A, et al. Resistance mechanisms to targeted therapies in ROS1(+) and ALK(+) non-small cell lung cancer. Clin Cancer Res 24, 3334-3347 (2018).
- 26. Lindeman, NI, Cagle, PT, Aisner, DL, Arcila, ME, Beasley, MB, Bernicker, EH, et al. Updated molecular testing guideline for the selection of lung cancer patients for treatment with targeted tyrosine kinase inhibitors: guideline from the College of American Pathologists, the International Association for the Study of Lung Cancer, and the Association for Molecular Pathology. J Thorac Oncol 13, 323-358 (2018).
- Ewalt, MD, West, H & Aisner, DL. Next generation sequencing-testing multiple genetic markers at once. JAMA Oncol 5, 1076 (2019).
- 28. Liam, CK, Mallawathantri, S & Fong, KM. Is tissue still the issue in detecting molecular alterations in lung cancer? *Respirology* **25**, 933-943 (2020).

#### **ACKNOWLEDGEMENTS**

The results published here are in part based upon data generated by the TCGA Research Network; https://www.cancer.gov/tcqa.

#### **AUTHOR CONTRIBUTIONS**

JZ and NPY conceived the project; CM and IB study design; CM, DEF, YM and RY data acquisition; CM, IH and NPY data analysis and interpretation; EO, IG and IB data interpretation; JZ software development; CM, EO, IG, NPY and IB writing and revising the manuscript.

#### **FUNDING**

Imagene Al.

#### **COMPETING INTERESTS**

IG, IH, JZ and NPY are employees of Imagene AI and have stock options. CM was reimbursed by Imagene AI for attending a conference. EO, DEF, YM, RY and IB declare no competing interests.

#### **ETHICS APPROVAL**

The study was approved by the Ethics Committee, IRB 7451-20-SMC, in accordance with the Declaration of Helsinki.

#### **ADDITIONAL INFORMATION**

**Supplementary information** The online version contains supplementary material available at https://doi.org/10.1038/s41379-022-01141-4.

Correspondence and requests for materials should be addressed to Chen Mayer.

Reprints and permission information is available at http://www.nature.com/reprints

**Publisher's note** Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

CC (i) Op

**Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing,

adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/.

© The Author(s) 2022